

## Anti-JAKMIP2 Antibody

Catalog Number: A13357-1

### About JAKMIP2

JAKMIP2, also known as neuroendocrine long coiled-coil protein 1 (NECC1), is a member of the JAKMIP family of proteins whose members are thought play a role in JAK1 signaling and the regulation of microtubule cytoskeleton rearrangements. JAKMIP2 is expressed predominantly in neuronal cells but unlike JAKMIP1, JAKMIP2 localizes to the Golgi complex. It has been suggested that JAKMIP2 may play important roles in the control of the regulated secretory pathway.

### Overview

Product Name	Anti-JAKMIP2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-JAKMIP2 Antibody (Catalog # A13357-1). Tested in ELISA, WB applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	JAKMIP2 Antibody is supplied in PBS containing 0.02% sodium azide.
Storage Instructions	JAKMIP2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Host	Rabbit
Uniprot ID	Q96AA8


### Technical Details

Immunogen	Rabbit polyclonal JAKMIP2 antibody was raised against a 16 amino acid peptide near the center of human JAKMIP2. The immunogen is located within amino acids 250 - 300 of JAKMIP2.
Predicted Reactive Species	Bovine
Cross Reactivity	At least four isoforms of JAKMIP2 are known to exist; this antibody will detect all four isoforms. JAKMIP2 antibody is predicted to not cross-react with JAKMIP1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/mL
Purification	JAKMIP2 Antibody is affinity chromatography purified via peptide column.
Suggested Dilutions	JAKMIP2 antibody can be used for detection of JAKMIP2 by Western blot at 1 - 2 ug/mL.

Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested. Optimal dilutions for each application should be determined by the researcher.

## Anti-JAKMIP2 Antibody (A13357-1) Images

---



Western blot analysis of JAKMIP2 in mouse brain tissue lysate with JAKMIP2 antibody at 1 ug/mL.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-JAKMIP2 Antibody

For Research Use Only. Not for use in diagnostic procedures.