

Anti-SH3BGRL Antibody

Catalog Number: A13363-1

About SH3BGRL

SH3BGRL, an adaptor protein, is one member of the SH3BGR family. SH3BGRL is ubiquitously expressed in varied human tissues and organs, including the bone marrow, heart, lung, liver, and kidney. SH3BGRL encodes a small protein with a conserved proline-rich PLPPQIF region and two motifs, Homer EVH1-binding and SH3-binding motifs. As a scaffold protein, SH3BGRL could participate in the protein-protein interaction for the integration or/and crosstalk of signal transduction, membrane trafficking, cytoskeletal rearrangements, and other key cellular processes.

Overview

Product Name	Anti-SH3BGRL Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-SH3BGRL Antibody catalog # A13363-1. Tested in WB, IHC, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O75368

Technical Details

Immunogen	E.coli-derived human SH3BGRL recombinant protein (Position: I3-K111).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SH3BGRL Antibody

For Research Use Only. Not for use in diagnostic procedures.