

Anti-PCTAIRE-2 CDK17 Antibody

Catalog Number: A13416


Overview

Product Name	Anti-PCTAIRE-2 CDK17 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PCTAIRE-2 CDK17 Antibody catalog # A13416. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q00537

Technical Details

Immunogen	Synthesized peptide derived from PCTAIRE-2 . at AA range: 70-150
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:20000

Anti-PCTAIRE-2 CDK17 Antibody (A13416) Images



Western blot analysis of PCTK2 Antibody. The lane on the right is blocked with the PCTK2 peptide.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PCTAIRE-2 CDK17 Antibody

For Research Use Only. Not for use in diagnostic procedures.