

Anti-Annexin A10 ANXA10 Antibody

Catalog Number: A13560


Overview

Product Name	Anti-Annexin A10 ANXA10 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Annexin A10 ANXA10 Antibody catalog # A13560. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UJ72

Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000

Anti-Annexin A10 ANXA10 Antibody (A13560) Images



Western blot analysis of Annexin A10 on different lysates using anti-Annexin A10 antibody at 1/1,000 dilution. Positive control: Lane 1: Hela Lane 2: Jurkat

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Annexin A10 ANXA10 Antibody

For Research Use Only. Not for use in diagnostic procedures.