

## Anti-PSP PSPN Antibody

Catalog Number: A13643

### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-PSP PSPN Antibody  |
| Reactive Species     | Human, Mouse, Rat   |
| Description          | Boster Bio Anti-PSP PSPN Antibody catalog # A13643. Tested in ELISA, WB applications. This antibody reacts with Human.  |
| Application          | ELISA, IHC, WB  |
| Clonality            | Polyclonal  |
| Formulation          | Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.  |
| Host                 | Rabbit  |
| Uniprot ID           | O60542  |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | The antiserum was produced against synthesized peptide derived from the C-terminal region of human PSPN. AA range:101-150 |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 1 mg/ml   |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.     |
| Suggested Dilutions | WB 1:500-2000<br>IHC 1:50-300<br>ELISA 1:2000-20000   |

## Anti-PSP PSPN Antibody (A13643) Images



Western blotting validation for Anti-PSP PSPN Antibody A13643 Western Blot (WB) analysis of 293 cells using PSP polyclonal antibody. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PSP PSPN Antibody

For Research Use Only. Not for use in diagnostic procedures.