

Anti-Centrosomal protein of 78 kDa CEP78 Antibody

Catalog Number: A13706-1


Overview

Product Name	Anti-Centrosomal protein of 78 kDa CEP78 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Centrosomal protein of 78 kDa CEP78 Antibody catalog # A13706-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q5JTW2

Technical Details

Immunogen	Synthetic peptide, corresponding Human CEP78.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Centrosomal protein of 78 kDa CEP78 Antibody (A13706-1) Images



Western blot (WB) analysis of CEP78 polyclonal antibody at 1:500 dilution Lane1: COS-7 whole cell lysate(40ug) Lane2: HEK293T whole cell lysate(40ug) Lane3: CT26 whole cell lysate(40ug) Lane4: H9C2 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Centrosomal protein of 78 kDa CEP78 Antibody

For Research Use Only. Not for use in diagnostic procedures.