

Anti-SRp46 (E72) SRSF8 Antibody

Catalog Number: A13713

Overview

| | |
|----------------------|--|
| Product Name | Anti-SRp46 (E72) SRSF8 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-SRp46 (E72) SRSF8 Antibody catalog # A13713. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat. |
| Application | IHC, WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q9BRL6 |

Technical Details

| | |
|---------------------|---|
| Immunogen | Synthetic peptide, corresponding to amino acids 50-100 of Human SRp46. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000 IHC: 1:50-1:200 |

Anti-SRp46 (E72) SRSF8 Antibody (A13713) Images



Western blot (WB) analysis of SRp46 (E72) pAb at 1:500 dilution Lane1:PC12 whole cell lysate(40ug) Lane2:The Embryo tissue lysate of Mouse(40ug) Lane3:EC9706 whole cell lysate(20ug) Lane4:HEK293T whole cell lysate(10ug) Lane5:A549 whole cell lysate(20ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SRp46 (E72) SRSF8 Antibody

For Research Use Only. Not for use in diagnostic procedures.