

Anti-BSX Antibody

Catalog Number: A13902

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-BSX Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-BSX Antibody catalog # A13902. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q3C1V8

Technical Details

Immunogen	Synthesized peptide derived from BSX . at AA range: 150-230
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:5000

Anti-BSX Antibody (A13902) Images



Western Blot analysis of 293T cells using BSX Polyclonal Antibody diluted at 1:500

Anti-BSX Antibody

For Research Use Only. Not for use in diagnostic procedures.