

## Anti-ADAM32 (N530) Antibody

Catalog Number: A13911-1


### Overview

Product Name	Anti-ADAM32 (N530) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ADAM32 (N530) Antibody catalog # A13911-1. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8TC27

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 500-560 of Human ADAM32.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-ADAM32 (N530) Antibody (A13911-1) Images



Western blot (WB) analysis of ADAM32 (N530) polyclonal antibody at 1:500 dilution Lane1:BV2 whole cell lysate(40ug) Lane2:The Testis tissue lysate of Rat(40ug) Lane3:MCF-7 whole cell lysate(40ug) Lane4:Hela whole cell lysate(20ug) Lane5:U-87MG whole cell lysate(30ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ADAM32 (N530) Antibody

For Research Use Only. Not for use in diagnostic procedures.