

## Anti-GPRC5D Antibody

Catalog Number: A13912


### Overview

Product Name	Anti-GPRC5D Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-GPRC5D Antibody catalog # A13912. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NZD1

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 280-330 of Human GPRC5D.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-GPRC5D Antibody (A13912) Images



Western blot (WB) analysis of GPRC5D pAb at 1:500 dilution  
Lane1:SGC7901 whole cell lysate(40ug) Lane2:PANC1 whole  
cell lysate(40ug) Lane3:CT26 whole cell lysate(40ug)  
Lane4:PC12 whole cell lysate(40ug)

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-GPRC5D Antibody

For Research Use Only. Not for use in diagnostic procedures.