

Anti-NIDD/ZDHHC23 Antibody

Catalog Number: A13972


Overview

Product Name	Anti-NIDD/ZDHHC23 Antibody
Reactive Species	Human
Description	Boster Bio Anti-NIDD/ZDHHC23 Antibody (Catalog# A13972). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	50 ul sera
Storage Instructions	Store the antibody at 4°C, stable for 6 months. For long-term storage, store at -20°C. Avoid repeated freeze and thaw cycles.
Host	Rabbit
Uniprot ID	Q8IYP9

Technical Details

Immunogen	A Full-length recombinant of human NIDD protein
Predicted Reactive Species	Mouse, Rat
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Suggested Dilutions	WB: 1:1000-1:4000

Anti-NIDD/ZDHHC23 Antibody (A13972) Images



Western blot analysis of NIDD. Anti-NIDD antibody was used with 1:4000 dilution on h Spleen lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NIDD/ZDHHC23 Antibody

For Research Use Only. Not for use in diagnostic procedures.