

Anti-K0746 Antibody (N-term)

Catalog Number: A13977

About SEL1L3

Orphan receptor.


Overview

Product Name	Anti-K0746 Antibody (N-term)
Reactive Species	Human
Description	Boster Bio Anti-K0746 Antibody (N-term) (Catalog # A13977). Tested in Flow Cytometry, WB application(s). This antibody reacts with Human.
Application	Flow Cytometry, WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q68CR1


Technical Details

Immunogen	This K0746 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 253-282 amino acids from the N-terminal region of human K0746.
Predicted Reactive Species	Mouse, Rat
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000 FC: 1:10-1:50

Anti-K0746 Antibody (N-term) (A13977) Images



All lanes : Anti-K0746 Antibody (N-term) at 1:1000 dilution
Lane 1: PC-3 whole cell lysate
Lane 2: U-251 MG whole cell lysate
Lane 3: U-87 MG whole cell lysate
Lane 4: A431 whole cell lysate
Lane 5: U-2OS whole cell lysate
Lysates/proteins at 20 µg per lane. Secondary Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/10000 dilution. Predicted band size : 125 kDa
Blocking/Dilution buffer: 5% NFDM/TBST.



K0746 Antibody (N-term) flow cytometric analysis of U251 cells (right histogram) compared to a negative control cell (left histogram). FITC-conjugated goat-anti-rabbit secondary antibodies were used for the analysis.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-K0746 Antibody (N-term)

For Research Use Only. Not for use in diagnostic procedures.