

## Anti-NEDD4-binding protein 2-like 2 N4BP2L2 Antibody

Catalog Number: A14016-1


### Overview

|                      |                                                                                                                                                           |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-NEDD4-binding protein 2-like 2 N4BP2L2 Antibody                                                                                                      |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                         |
| Description          | Boster Bio Anti-NEDD4-binding protein 2-like 2 N4BP2L2 Antibody catalog # A14016-1. Tested in WB applications. This antibody reacts with Human,Mouse,Rat. |
| Application          | WB                                                                                                                                                        |
| Clonality            | Polyclonal                                                                                                                                                |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2                                                                                    |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.                |
| Host                 | Rabbit                                                                                                                                                    |
| Uniprot ID           | Q92802                                                                                                                                                    |

### Technical Details

|                     |                                                                                                                                                             |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen           | Synthetic peptide, corresponding to amino acids 400-445 of Human N4BP2L2.                                                                                   |
| Isotype             | IgG                                                                                                                                                         |
| Form                | Liquid                                                                                                                                                      |
| Concentration       | 1 mg/ml                                                                                                                                                     |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000                                                                                                                                            |

## Anti-NEDD4-binding protein 2-like 2 N4BP2L2 Antibody (A14016-1) Images



Western blot (WB) analysis of N4BP2L2 polyclonal antibody at 1:500 dilution Lane1:MCF-7 cell lysate Lane2:H9C2 cell lysate Lane3:Rat lung tissue lysate

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NEDD4-binding protein 2-like 2 N4BP2L2 Antibody

For Research Use Only. Not for use in diagnostic procedures.