

Anti-FAM80B RIMKLB Antibody

Catalog Number: A14138


Overview

Product Name	Anti-FAM80B RIMKLB Antibody
Reactive Species	Human, Monkey, Mouse, Rat
Description	Boster Bio Anti-FAM80B RIMKLB Antibody catalog # A14138. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat, Monkey.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9ULI2


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human RIMKB. AA range:177-226
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:5000

Anti-FAM80B RIMKLB Antibody (A14138) Images



Western Blot analysis of various cells using FAM80B Polyclonal Antibody diluted at 1:500



Western blot analysis of lysates from HepG2 and COS cells, using RIMKB Antibody. The lane on the right is blocked with the synthesized peptide.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FAM80B RIMKLB Antibody

For Research Use Only. Not for use in diagnostic procedures.