

## Anti-WDR32/DCAF10 Antibody

Catalog Number: A14141-1

### About DCAF10

WDR32, containing seven WD repeats, is also termed DDB1- and CUL4-associated factor 10 (DCAF10), which may function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex (PMID: 16949367). WDR32 is able to interact with DDB1 (PMID: 16949367). WDR32 is supposed to be implicated in process of ubiquitin-dependent protein degradation system. Two isoforms are described and produced by alternative splicing. Proteintech's WDR32 antibody 20483-1-AP majorly detects WDR32 isoform1 with a theoretical MW of 61 kDa.

### Overview

Product Name	Anti-WDR32/DCAF10 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-WDR32/DCAF10 Antibody catalog # A14141-1. Tested in WB, IHC, ICC, IF, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q5QP82

### Technical Details

Immunogen	E.coli-derived human WDR32/DCAF10 recombinant protein (Position: F3-P223).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ImmunoPrecipitation, 1:250-300

ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-WDR32/DCAF10 Antibody

For Research Use Only. Not for use in diagnostic procedures.