

## Anti-Cytochrome P450 20A1 CYP20A1 Antibody

Catalog Number: A14222-1

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




### Overview

|                      |                                                                                                                                                   |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-Cytochrome P450 20A1 CYP20A1 Antibody                                                                                                        |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                 |
| Description          | Boster Bio Anti-Cytochrome P450 20A1 CYP20A1 Antibody catalog # A14222-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat. |
| Application          | WB                                                                                                                                                |
| Clonality            | Polyclonal                                                                                                                                        |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2                                                                            |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.        |
| Host                 | Rabbit                                                                                                                                            |
| Uniprot ID           | Q6UW02                                                                                                                                            |

### Technical Details

|                     |                                                                                                                                                             |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen           | Synthetic peptide, corresponding to amino acids 380-430 of Human CYP20A1.                                                                                   |
| Isotype             | IgG                                                                                                                                                         |
| Form                | Liquid                                                                                                                                                      |
| Concentration       | 1 mg/ml                                                                                                                                                     |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000                                                                                                                                            |

## Anti-Cytochrome P450 20A1 CYP20A1 Antibody (A14222-1) Images



Western blot (WB) analysis of CYP20A1 polyclonal antibody at 1:500 dilution Lane1:MCF-7 whole cell lysate(40ug) Lane2:HepG2 whole cell lysate(40ug) Lane3:AML-12 whole cell lysate(40ug) Lane4:PC12 whole cell lysate(40ug)

[Anti-Cytochrome P450 20A1 CYP20A1 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.