

Anti-NAP1L5 Antibody

Catalog Number: A14492

About NAP1L5

The nucleosome assembly protein 1-like 5 (NAP1L5) protein, also known as DRLM, was initially identified as a possible tumor suppressor gene in hepatocellular carcinoma (1). The NAP1L5 gene is located within the intron of Herc3, and unlike Herc3, is imprinted and highly expressed in regions of the adult mouse brain (2).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-NAP1L5 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NAP1L5 Antibody (Catalog # A14492). Tested in ELISA, WB, IHC-P, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC-P, WB
Clonality	Polyclonal
Formulation	NAP1L5 antibody is supplied in PBS containing $\leq 0.09\%$ sodium azide.
Storage Instructions	NAP1L5 antibody can be stored at 4°C, stable for up to one year.
Host	Rabbit
Uniprot ID	Q96NT1

Technical Details


Immunogen	NAP1L5 antibody was raised against an 18 amino acid peptide near the center of human NAP1L5. The immunogen is located within amino acids 40 - 90 of NAP1L5.
Predicted Reactive Species	Bovine
Cross Reactivity	NAP1L5 antibody is human, mouse and rat reactive.
Isotype	IgG

Form	Liquid
Concentration	1 mg/mL
Purification	NAP1L5 antibody is affinity chromatography purified via peptide column.
Suggested Dilutions	NAP1L5 antibody can be used for detection of NAP1L5 by Western blot at 1 - 2 ug/mL. Antibody can also be used for immunohistochemistry starting at 5 ug/mL. For immunofluorescence start at 20 ug/mL. Antibody validated: Western Blot in mouse samples; Immunohistochemistry in mouse samples and Immunofluorescence in mouse samples. All other applications and species not yet tested. Optimal dilutions for each application should be determined by the researcher.


Anti-NAP1L5 Antibody (A14492) Images



Western blot analysis of NAP1L5 in EL4 cell lysate with NAP1L5 antibody at 1 ug/ml in (A) the absence and (B) the presence of blocking peptide.



Immunohistochemistry of NAP1L5 in mouse brain tissue with NAP1L5 antibody at 5 ug/mL.



Immunofluorescence of NAP1L5 in mouse brain tissue with NAP1L5 antibody at 20 ug/mL.

Anti-NAP1L5 Antibody

For Research Use Only. Not for use in diagnostic procedures.