

Anti-Histone H4 HIST1H4A Antibody

Catalog Number: A14495-2


Overview

Product Name	Anti-Histone H4 HIST1H4A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Histone H4 HIST1H4A Antibody catalog # A14495-2. Tested in WB,ICC/IF,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P62805


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 ICC: 1:50-1:200 IHC: 1:50-1:200

Anti-Histone H4 HIST1H4A Antibody (A14495-2) Images



Western blot analysis of Histone H4 on mouse spleen lysates using anti-Histone H4 antibody at 1/5,000 dilution.



ICC staining Histone H4 in PANC-1 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Histone H4 HIST1H4A Antibody

For Research Use Only. Not for use in diagnostic procedures.