

## Anti-ACSS3 Antibody

Catalog Number: A14541

### About ACSS3

May be involved in transcriptional regulation.

Lu X.W., Zhongguo Yi Xue Ke Xue Yuan Xue Bao 23:27-31(2001).

Xu X., Submitted (MAY-2000) to the EMBL/GenBank/DDBJ databases.

Ota T., Nat. Genet. 36:40-45(2004).

### Overview

Product Name	Anti-ACSS3 Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-ACSS3 Antibody catalog # A14541. Tested in WB applications. This antibody reacts with Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H6R3

### Technical Details

Immunogen	Synthesized peptide derived from internal of human TCEAL4.
Predicted Reactive Species	Chimpanzee, Drosophila, Macaque
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.

If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.  
Some PubMed article(s) citing the expression level of this target are as follows:  
Boster Bio's internal QC testing used:  
WB: 1:500-1:2000

## Anti-ACSS3 Antibody (A14541) Images




Figure 1. Western blotting validation for Anti-ACSS3 Antibody A14541

Western blot analysis of extracts of various cell lines

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ACSS3 Antibody