

## Anti-TP8L1 TNFAIP8L1 Antibody

Catalog Number: A14573-1


### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-TP8L1 TNFAIP8L1 Antibody  |
| Reactive Species     | Human, Mouse   |
| Description          | Boster Bio Anti-TP8L1 TNFAIP8L1 Antibody catalog # A14573-1. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse.     |
| Application          | ELISA, WB  |
| Clonality            | Polyclonal   |
| Formulation          | Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.  |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit   |
| Uniprot ID           | Q8WVP5   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | Synthesized peptide derived from human protein at AA range: 40-120  |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 1 mg/ml   |
| Purification        | TP8L1 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Suggested Dilutions | WB 1:500-2000<br>ELISA 1:5000-20000   |

## Anti-TP8L1 TNFAIP8L1 Antibody (A14573-1) Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4°over night

## 2 Publications Citing This Product

1. PubMed ID: -, Zhang X, Li Z, Sun Y, Liu G, Liu X, Zhou Y. Overexpression of Tumour Necrosis Factor-alpha-Induced Protein 8 is Associated with Prognosis in Colon Cancer. *Cancer Manag Res.* 2021;13:4055-4065 <https://doi.org/10.2147/CMAR.S297451>
2. PubMed ID: 31929755, Mengya Zhong, Nini Li, Xingfeng Qiu, Yuhan Ye, Huiyu Chen, Jianyu Hua, Ping Yin, Guohong Zhuang *Int J Biol Sci.* 2020; 16(2): 272-283. Published online 2020 Jan 1. doi:10.7150/ijbs.37906

Visit [bosterbio.com/anti-tp8l1-tnfaip8l1-antibody-a14573-1-boster.html](https://www.bosterbio.com/anti-tp8l1-tnfaip8l1-antibody-a14573-1-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TP8L1 TNFAIP8L1 Antibody

For Research Use Only. Not for use in diagnostic procedures.