

## Anti-ARMCX3 Antibody

Catalog Number: A14662-2

### About ARMCX3

ARMCX3, also named as Armadillo repeat-containing X-linked protein 3, is a 379 amino acid protein, which contains 3 ARM repeats. ARMCX3 is a single pass membrane protein belonging to the armadillo repeat family of proteins. ARMCX3 is believed to play a role in embryonic development and tissue maintenance and may also function as a tumor suppressor.

### Overview

Product Name	Anti-ARMCX3 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-ARMCX3 Antibody catalog # A14662-2. Tested in WB, IHC, IP, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IP, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9UH62

### Technical Details

Immunogen	E.coli-derived human ARMCX3 recombinant protein (Position: E40-L143).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-ARMCX3 Antibody

For Research Use Only. Not for use in diagnostic procedures.