

## Anti-FAM92A1/CIBAR1 Antibody

Catalog Number: A14866

### About CIBAR1

Enables phospholipid binding activity. Involved in several processes, including inner mitochondrial membrane organization; limb morphogenesis; and membrane tubulation. Located in several cellular components, including centriole; ciliary base; and mitochondrial crista.

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-FAM92A1/CIBAR1 Antibody                                                                                                                                                                                                                                                                                                                                                                                               |
| Reactive Species     | Human                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Description          | Boster Bio Anti-FAM92A1/CIBAR1 Antibody catalog # A14866. Tested in WB, IHC, ICC, IF, IP, ELISA applications. This antibody reacts with Human.                                                                                                                                                                                                                                                                             |
| Application          | ELISA, IP, IF, IHC, ICC, WB                                                                                                                                                                                                                                                                                                                                                                                                |
| Clonality            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Formulation          | 500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol<br>This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used.<br>For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | 12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.                                                                                                                                                                                                                                                                                          |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Uniprot ID           | A1XB55                                                                                                                                                                                                                                                                                                                                                                                                                     |

### Technical Details

|               |                                                                                |
|---------------|--------------------------------------------------------------------------------|
| Immunogen     | E.coli-derived human FAM92A1/CIBAR1 recombinant protein (Position: L164-T316). |
| Form          | Liquid                                                                         |
| Concentration | 500 ug/ml                                                                      |

|                     |                                                                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Purification        | Immunogen affinity purified.                                                                                                                                          |
| Suggested Dilutions | Western blot, 1:500-2000<br>Immunohistochemistry, 1:50-400<br>Immunocytochemistry/Immunofluorescence, 1:50-400<br>ImmunoPrecipitation, 1:250-300<br>ELISA, 1:100-1000 |

[Anti-FAM92A1/CIBAR1 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.