

Anti-Protein ABHD12B ABHD12B Antibody

Catalog Number: A14869


Overview

Product Name	Anti-Protein ABHD12B ABHD12B Antibody
Reactive Species	Human
Description	Boster Bio Anti-Protein ABHD12B ABHD12B Antibody catalog # A14869. Tested in IHC applications. This antibody reacts with Human.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q7Z5M8

Technical Details

Immunogen	A peptide derived from human ABHD12B. Immunogen sequence location: 281-330
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used User need to optimize the dilution ratio for this antibody For protocols please visit https://www.bosterbio.com/protocol-and-troubleshooting/

Anti-Protein ABHD12B ABHD12B Antibody (A14869) Images



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma, using ABHD12B Antibody. The picture on the right is blocked with the synthesized peptide.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Protein ABHD12B ABHD12B Antibody