

## Anti-HSP77/76 Antibody

Catalog Number: A14896

### Overview


Product Name	Anti-HSP77/76 Antibody
Reactive Species	Human
Description	Boster Bio Anti-HSP77/76 Antibody catalog # A14896. Tested in WB, ELISA applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P48741

### Technical Details

Immunogen	Synthesized peptide derived from the Internal region of human HSP77/76.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:20000

## Anti-HSP77/76 Antibody (A14896) Images

---



Western Blot analysis of K562 cells using HSP77/76 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-HSP77/76 Antibody

For Research Use Only. Not for use in diagnostic procedures.