

Anti-ALS2CR13 FAM117B Antibody

Catalog Number: A14981


Overview

| | |
|----------------------|---|
| Product Name | Anti-ALS2CR13 FAM117B Antibody |
| Reactive Species | Human, Mouse |
| Description | Boster Bio Anti-ALS2CR13 FAM117B Antibody catalog # A14981. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse. |
| Application | ELISA, WB |
| Clonality | Polyclonal |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q6P1L5 |


Technical Details

| | |
|---------------------|--|
| Immunogen | The antiserum was produced against synthesized peptide derived from human ALS2CR13. AA range:411-460 |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | Immunogen affinity purified |
| Suggested Dilutions | WB 1:500-1:2000 ELISA 1:5000 |

Anti-ALS2CR13 FAM117B Antibody (A14981) Images



Western Blot analysis of HELA cells using ALS2CR13 Polyclonal Antibody diluted at 1:500



Western Blot analysis of HELA cells using ALS2CR13 Polyclonal Antibody diluted at 1:500

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ALS2CR13 FAM117B Antibody

For Research Use Only. Not for use in diagnostic procedures.