

## Anti-RPL22L1 Antibody Picoband®

Catalog Number: A15161-2

### About RPL22L1

Predicted to enable RNA binding activity. Predicted to be a structural constituent of ribosome. Predicted to be involved in cytoplasmic translation. Predicted to be located in ribosome. Predicted to be part of ribonucleoprotein complex.

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview


|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-RPL22L1 Antibody Picoband®                                                                                                                                                                                                                                                                                                                                                                                                          |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Description          | Boster Bio Anti-RPL22L1 Antibody Picoband® catalog # A15161-2. Tested in WB, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance. |
| Application          | ELISA, IF, ICC, WB                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Clonality            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Formulation          | Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na <sub>2</sub> HPO <sub>4</sub> .                                                                                                                                                                                                                                                                                                                                                |
| Storage Instructions | At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.                                                                                                                                                                                                                                        |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Uniprot ID           | Q6P5R6                                                                                                                                                                                                                                                                                                                                                                                                                                   |

### Technical Details


|               |                                                                           |
|---------------|---------------------------------------------------------------------------|
| Immunogen     | E.coli-derived human RPL22L1 recombinant protein (Position: M1-D122).     |
| Form          | Lyophilized                                                               |
| Concentration | Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml. |

|                     |                                                                                                                                   |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Purification        | Immunogen affinity purified.                                                                                                      |
| Suggested Dilutions | Western blot, 0.25-0.5 ug/ml, Human, Mouse, Rat<br>Immunocytochemistry/Immunofluorescence, 5 ug/ml, Human<br>ELISA, 0.1-0.5 ug/ml |

## Anti-RPL22L1 Antibody Picoband® (A15161-2) Images



Western blot analysis of RPL22L1 using anti-RPL22L1 antibody (A15161-2). Electrophoresis was performed on a 12% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 2 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human RT4 whole cell lysates, Lane 2: human A549 whole cell lysates, Lane 3: rat RH35 whole cell lysates, Lane 4: rat C6 whole cell lysates, Lane 5: mouse HEPA1-6 whole cell lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-RPL22L1 antigen affinity purified polyclonal antibody (A15161-2) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196-200) with Tanon 5200 system. A specific band was detected for RPL22L1 at approximately 15 kDa. The expected band size for RPL22L1 is at 15 kDa.



IF analysis of RPL22L1 using anti-RPL22L1 antibody (A15161-2). RPL22L1 was detected in an immunocytochemical section of A549 cells. Enzyme antigen retrieval was performed using IHC enzyme antigen retrieval reagent (AR0022) for 15 mins. The cells were blocked with 10% goat serum. And then incubated with 5 ug/mL rabbit anti-RPL22L1 Antibody (A15161-2) overnight at 4°C. Cy3 Conjugated Goat Anti-Rabbit IgG (BA1032) was used as secondary antibody at 1:500 dilution and incubated for 30 minutes at 37°C. The section was counterstained with DAPI. Visualize using a fluorescence microscope and filter sets appropriate for the label used.

### Anti-RPL22L1 Antibody

For Research Use Only. Not for use in diagnostic procedures.