

## Anti-UBQLNL Antibody

Catalog Number: A15200

### About UBQLNL

UBQLNL also named as Ubiquilin like protein, encodes a 475 amino acid protein. The specific function of it is still non-known. UBQLNL may have an association with type 2 diabetes, but the mechanism is supposed be studied in the further research. The subcellular localization of UBQLNL is nucleus. The expression of it is detected in some tissues, such as brain ,testis, heart and liver.

### Overview

Product Name	Anti-UBQLNL Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-UBQLNL Antibody catalog # A15200. Tested in WB, IHC, ICC, IF, IP, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8IYU4

### Technical Details

Immunogen	E.coli-derived human UBQLNL recombinant protein (Position: R134-L473).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-UBQLNL Antibody

For Research Use Only. Not for use in diagnostic procedures.