

## Anti-TMCC2 Antibody

Catalog Number: A15276-1

### About TMCC2

TMCC2 also termed as Cerebral protein 11 or HUCEP11 is a 709 amino acid protein, which is a multi-pass membrane protein and belongs to the TEX28 family. The TMCC2 may have an association with Alzheimer's disease.

### Overview

Product Name	Anti-TMCC2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TMCC2 Antibody catalog # A15276-1. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O75069

### Technical Details

Immunogen	E.coli-derived human TMCC2 recombinant protein (Position: R3-L230).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TMCC2 Antibody

For Research Use Only. Not for use in diagnostic procedures.