

Anti-Olfactory receptor 2A25 Antibody

Catalog Number: A15315


Overview

| | |
|----------------------|---|
| Product Name | Anti-Olfactory receptor 2A25 Antibody |
| Reactive Species | Human |
| Description | Boster Bio Anti-Olfactory receptor 2A25 Antibody catalog # A15315. Tested in WB, IF, ELISA applications. This antibody reacts with Human. |
| Application | ELISA, IF, WB |
| Clonality | Polyclonal |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | A4D2G3 |


Technical Details

| | |
|---------------------|--|
| Immunogen | The antiserum was produced against synthesized peptide derived from human OR2A25. AA range:241-290 |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | Immunogen affinity purified |
| Suggested Dilutions | WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:5000 |


Anti-Olfactory receptor 2A25 Antibody (A15315) Images




Western Blot analysis of various cells using Olfactory receptor 2A25 Polyclonal Antibody



Western blot analysis of lysates from HT-29 cells, using OR2A25 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using OR2A25 antibody.



Immunofluorescence analysis of A549 cells, using OR2A25 Antibody. The picture on the right is blocked with the synthesized peptide.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Olfactory receptor 2A25 Antibody

For Research Use Only. Not for use in diagnostic procedures.