

Anti-SUHW1 ZNF280A Antibody

Catalog Number: A15519-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-SUHW1 ZNF280A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-SUHW1 ZNF280A Antibody catalog # A15519-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P59817

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 491-536 of Human SUHW1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-SUHW1 ZNF280A Antibody (A15519-1) Images



Western blot (WB) analysis of SUHW1 polyclonal antibody at 1:500 dilution Lane1:HEK293T cell lysate Lane2:NIH-3T3 cell lysate Lane3:Rat kidney tissue lysate

Anti-SUHW1 ZNF280A Antibody

For Research Use Only. Not for use in diagnostic procedures.