

Anti-R Cbwd1 Antibody (C-term)

Catalog Number: A15555

About ZNG1

This gene encodes a member of the type I (acidic) keratin family, which belongs to the superfamily of intermediate filament (IF) proteins. Keratins are heteropolymeric structural proteins which form the intermediate filament. These filaments, along with actin microfilaments and microtubules, compose the cytoskeleton of epithelial cells. The type I keratin genes are clustered in a region of chromosome 17q12-q21.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-R Cbwd1 Antibody (C-term)
Reactive Species	Rat
Description	Boster Bio Anti-R Cbwd1 Antibody (C-term) (Catalog # A15555). Tested in WB application(s). This antibody reacts with Rat.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99MB4


Technical Details

Immunogen	This Rat Cbwd1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 325-351 amino acids from the C-terminal region of rat Cbwd1.
Predicted Reactive Species	Mouse, Rat
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.

Suggested Dilutions

WB: 1:1000

Anti-R Cbwd1 Antibody (C-term) (A15555) Images



Rat Cbwd1 Antibody (C-term) western blot analysis in rat bladder, cerebellum and stomach tissue lysates (35ug/lane). This demonstrates the Rat Cbwd1 antibody detected the Rat Cbwd1 protein (arrow).

Anti-R Cbwd1 Antibody (C-term)

For Research Use Only. Not for use in diagnostic procedures.