

Anti-TEX264 Antibody

Catalog Number: A15581-2

About TEX264

Enables signaling receptor activity. Involved in protein-DNA covalent cross-linking repair. Acts upstream of or within reticulophagy. Located in several cellular components, including autophagosome membrane; endoplasmic reticulum membrane; and replication fork.

Overview

Product Name	Anti-TEX264 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-TEX264 Antibody catalog # A15581-2. Tested in WB, IHC, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9Y6I9

Technical Details

Immunogen	E.coli-derived human TEX264 recombinant protein (Position: 36-310).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TEX264 Antibody

For Research Use Only. Not for use in diagnostic procedures.