

## Anti-C1orf26/SWT1 Antibody Picoband® FITC Conjugated

Catalog Number: A15590-FITC

### About SWT1

Predicted to be active in nucleus.

### Overview

Product Name	Anti-C1orf26/SWT1 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5T5J6

### Technical Details

Immunogen	E.coli-derived human C1orf26/SWT1 recombinant protein (Position: Q322-R889). Human C1orf26/SWT1 shares 80.7% amino acid (aa) sequence identity with mouse C1orf26/SWT1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C1orf26/SWT1 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.