

Anti-C1orf26/SWT1 Antibody Picoband® HRP Conjugated

Catalog Number: A15590-HRP

About SWT1

Predicted to be active in nucleus.

Overview

Product Name	Anti-C1orf26/SWT1 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q5T5J6

Technical Details

Immunogen	E.coli-derived human C1orf26/SWT1 recombinant protein (Position: Q322-R889). Human C1orf26/SWT1 shares 80.7% amino acid (aa) sequence identity with mouse C1orf26/SWT1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-C1orf26/SWT1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.