

Anti-C1orf26/SWT1 Antibody Picoband® PE Conjugated

Catalog Number: A15590-PE

About SWT1

Predicted to be active in nucleus.

Overview

Product Name	Anti-C1orf26/SWT1 Antibody Picoband® PE Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5T5J6

Technical Details

Immunogen	E.coli-derived human C1orf26/SWT1 recombinant protein (Position: Q322-R889). Human C1orf26/SWT1 shares 80.7% amino acid (aa) sequence identity with mouse C1orf26/SWT1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C1orf26/SWT1 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.