

Anti-SLC22A10 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A15610-1-Fluoro550

About SLC22A10

Solute carrier family 22 member 10 (SLC22A10), also known as organic anion transporter 5 (OAT5), is a protein that in humans is encoded by the SLC22A10 gene.

Overview

Product Name	Anti-SLC22A10 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q63ZE4

Technical Details

Immunogen	E.coli-derived human SLC22A10 recombinant protein (Position: N290-A541).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SLC22A10 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.