

Anti-HDHD2 Antibody Picoband® Biotin Conjugated

Catalog Number: A15632-1-Biotin

About HDHD2

Haloacid dehalogenase-like hydrolase domain-containing protein 2 is an enzyme that in humans is encoded by the HDHD2 gene. HDHD2 (haloacid dehalogenase-like hydrolase domain containing 2) is also known as DKFZp564D1378 and is a 259 amino acid protein that is expressed as two isoforms produced by alternative splicing. HDHD2 belongs to the HAD-like hydrolase superfamily, which contains a group of hydrolase enzymes.

Overview

Product Name	Anti-HDHD2 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9H0R4

Technical Details

Immunogen	E.coli-derived human HDHD2 recombinant protein (Position: M1-L259). Human HDHD2 shares 89.2% and 90% amino acid (aa) sequence identity with mouse and rat HDHD2, respectively.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HDHD2 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.