

## Anti-PERM1 Antibody Picoband® FITC Conjugated

Catalog Number: A15642-FITC

### About PERM1

PERM1 is predicted to play a role in glucose and lipid metabolism, energy transfer, and muscle contraction. In mouse muscle cells, Perm1 functions downstream of peroxisome proliferator-activated receptor-gamma coactivator-1 and estrogen-related receptors and regulates a subset of Pgc1- and Err-dependent genes.

### Overview

Product Name	Anti-PERM1 Antibody Picoband® FITC Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5SV97

### Technical Details

Immunogen	E.coli-derived human PERM1 recombinant protein (Position: R69-N763). Human PERM1 shares 50.1% amino acid (aa) sequence identity with mouse PERM1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-PERM1 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.