

Anti-GRAMD1A Antibody Picoband® FITC Conjugated

Catalog Number: A15650-FITC

About GRAMD1A

Predicted to enable cholesterol binding activity and cholesterol transfer activity. Predicted to be involved in cellular response to cholesterol and intracellular sterol transport. Located in several cellular components, including cytosol; endoplasmic reticulum membrane; and organelle membrane contact site.

Overview

Product Name	Anti-GRAMD1A Antibody Picoband® FITC Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q96CP6

Technical Details

Immunogen	E.coli-derived human GRAMD1A recombinant protein (Position: E293-S724).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GRAMD1A Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.