

## Anti-Zinc finger protein 597 ZNF597 Antibody

Catalog Number: A15681


### Overview

|                      |                                                                                                                                                        |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-Zinc finger protein 597 ZNF597 Antibody                                                                                                           |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                      |
| Description          | Boster Bio Anti-Zinc finger protein 597 ZNF597 Antibody catalog # A15681. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat. |
| Application          | IF, IHC, WB                                                                                                                                            |
| Clonality            | Polyclonal                                                                                                                                             |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2                                                                                 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.             |
| Host                 | Rabbit                                                                                                                                                 |
| Uniprot ID           | Q96LX8                                                                                                                                                 |

### Technical Details

|                     |                                                                                                                                                             |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen           | Recombinant protein of human ZNF597                                                                                                                         |
| Isotype             | IgG                                                                                                                                                         |
| Form                | Liquid                                                                                                                                                      |
| Concentration       | 1 mg/ml                                                                                                                                                     |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:2000<br>IHC: 1:50-1:200<br>IF: 1:50-1:200                                                                                                       |

## Anti-Zinc finger protein 597 ZNF597 Antibody (A15681) Images



Western blot analysis of extracts of various cell lines, using ZNF597 antibody.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zinc finger protein 597 ZNF597 Antibody

For Research Use Only. Not for use in diagnostic procedures.