

## Anti-SLC25A44 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A15684-1-Fluoro550

### About SLC25A44

SLC25A44 belongs to the SLC25 family of mitochondrial carrier proteins (Haitina et al., 2006 [PubMed 16949250]).

### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-SLC25A44 Antibody Picoband® Fluoro550 Conjugated   |
| Reactive Species     | Human, Mouse, Rat   |
| Application          | Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality            | Polyclonal  |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.  |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.  |
| Host                 | Rabbit  |
| Uniprot ID           | Q96H78  |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | A synthetic peptide corresponding to a sequence in the middle region of human SLC25A44. |
| Form                | Liquid  |
| Concentration       | 0.5 mg/mL   |
| Purification        | Immunogen affinity purified.  |
| Conjugate           | Fluoro550<br>Excitation Wavelength: 562 nm<br>Emission Wavelength: 576 nm               |
| Suggested Dilutions | Optimal dilutions should be determined by end users.                                    |

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SLC25A44 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.