

Anti-DCAF16 Antibody Picoband® Biotin Conjugated

Catalog Number: A15723-1-Biotin

About DCAF16

Predicted to be involved in protein ubiquitination. Part of Cul4-RING E3 ubiquitin ligase complex.

Overview

Product Name	Anti-DCAF16 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9NXF7

Technical Details

Immunogen	E.coli-derived human DCAF16 recombinant protein (Position: M1-L216).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.