

## Anti-DCAF16 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A15723-1-Fluoro488

### About DCAF16

Predicted to be involved in protein ubiquitination. Part of Cul4-RING E3 ubiquitin ligase complex.

### Overview

Product Name	Anti-DCAF16 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NXF7

### Technical Details

Immunogen	E.coli-derived human DCAF16 recombinant protein (Position: M1-L216).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DCAF16 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.