

Anti-DCAF16 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A15723-1-Fluoro594

About DCAF16

Predicted to be involved in protein ubiquitination. Part of Cul4-RING E3 ubiquitin ligase complex.

Overview

Product Name	Anti-DCAF16 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NXF7

Technical Details

Immunogen	E.coli-derived human DCAF16 recombinant protein (Position: M1-L216).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DCAF16 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.