

## Anti-DCAF16 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A15723-1-Fluoro594

### About DCAF16

Predicted to be involved in protein ubiquitination. Part of Cul4-RING E3 ubiquitin ligase complex.

### Overview

Product Name	Anti-DCAF16 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NXF7

### Technical Details

Immunogen	E.coli-derived human DCAF16 recombinant protein (Position: M1-L216).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DCAF16 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.