

## Anti-PRPF39 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A15727-2-Fluoro647

### About PRPF39

Predicted to be involved in mRNA 5'-splice site recognition. Predicted to be part of U1 snRNP; U2-type prespliceosome; and commitment complex.

### Overview

Product Name	Anti-PRPF39 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q86UA1

### Technical Details

Immunogen	E.coli-derived human PRPF39 recombinant protein (Position: D149-Y664).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-PRPF39 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.