

Anti-Olfactory receptor 4A15 Antibody

Catalog Number: A15816

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Olfactory receptor 4A15 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Olfactory receptor 4A15 Antibody catalog # A15816. Tested in WB, IF, ELISA applications. This antibody reacts with Human.
Application	ELISA, IF, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NGL6

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human OR4A15. AA range:261-310
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified


Suggested Dilutions

WB 1:500-1:2000
IF 1:200-1:1000
ELISA 1:10000


Anti-Olfactory receptor 4A15 Antibody (A15816) Images




Western Blot analysis of various cells using Olfactory receptor 4A15 Polyclonal Antibody



Western blot analysis of lysates from COLO and HepG2 cells, using OR4A15 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using OR4A15 antibody.



Immunofluorescence analysis of A549 cells, using OR4A15 Antibody. The picture on the right is blocked with the synthesized peptide.

[Anti-Olfactory receptor 4A15 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.