

## Anti-NTMT2 Antibody Picoband® Cy3 Conjugated

Catalog Number: A15877-1-Cy3

### About NTMT2

Enables N-terminal protein N-methyltransferase activity. Involved in N-terminal protein amino acid methylation. Located in nucleus.

### Overview

Product Name	Anti-NTMT2 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5VVY1

### Technical Details

Immunogen	E.coli-derived human NTMT2 recombinant protein (Position: M1-S283). Human NTMT2 shares 93.6% and 93.3% amino acid (aa) sequence identity with mouse and rat NTMT2, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NTMT2 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.