

## Anti-NTMT2 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A15877-1-Fluoro550

### About NTMT2

Enables N-terminal protein N-methyltransferase activity. Involved in N-terminal protein amino acid methylation. Located in nucleus.

### Overview

Product Name	Anti-NTMT2 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5VVY1

### Technical Details

Immunogen	E.coli-derived human NTMT2 recombinant protein (Position: M1-S283). Human NTMT2 shares 93.6% and 93.3% amino acid (aa) sequence identity with mouse and rat NTMT2, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NTMT2 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.