

Anti-AKAP14 Antibody

Catalog Number: A15886-3

Overview

| | |
|----------------------|--|
| Product Name | Anti-AKAP14 Antibody |
| Reactive Species | Human |
| Description | Boster Bio Anti-AKAP14 Antibody catalog # A15886-3. Tested in WB, IHC, ELISA applications. This antibody reacts with Human. |
| Application | ELISA, IHC, WB |
| Clonality | Polyclonal |
| Formulation | 500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | 12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | Q86UN6 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human AKAP14 recombinant protein (Position: M1-N195). |
| Form | Liquid |
| Concentration | 500 ug/ml |
| Purification | Immunogen affinity purified. |
| Suggested Dilutions | Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ELISA, 1:100-1000 |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-AKAP14 Antibody

For Research Use Only. Not for use in diagnostic procedures.